

### 3 MEMORY MODES

Enjoy your personal massage program according to your individual needs and expectations. The next time you use the chair, press the memory button and activate the automatic memo mode to enjoy your personal massage.

### MULTIPLE POSITIONS MASSAGE

Innovative acupressure point's sensor allows you to enjoy body massage in all parts of your body. By combining a variety of innovative methods, you can relax and reduce stress.

### 5 LEVELS OF MASSAGE WIDTH ADJUSTABLE

According to different body type, the back massage curve and shoulder position are automatically detected, and the massage range is adjusted according to different parts of the human body. When massaged under tapping and finger pressure, the massage width can be adjusted freely in the range of 6-21 CM with 5 levels.

# 5 LEVELS OF MASSAGE SPEED ADJUSTABLE

According to different massage needs, the massage speed can be selected in 5 levels. You can freely adjust the massage speed, design your massage according to your wishes and combine different massage experiences with exquisite quality to feel completely at ease.

### 3-LEVEL AIR PRESSURE INTENSITY ADJUSTABLE

With a total of 38 airbags and air pressure intensity adjustable in three levels, full relaxation from the shoulders to the feet

### LED TOUCHSCREEN

7 inches super touch screen, equipped with new UI operate system, making the massage progress clear to be seen. Users can control the massage program by one button.

# SKYLINER II – DRAWS YOU INTO A NEVER-BEFORE-SEEN WELLNESS SENSATION.



# TECHNICAL DATA

Dimensions upright $80 \times 148 \times 120 \text{ cm}$ Dimensions reclining $80 \times 180 \times 105 \text{ cm}$ 

Weight 141 kg

**Voltage** 220-240V ~ 50/60Hz

Nominal power 220 W
Rated runtime 20 minutes

braintronics®

Certificates



21 minutes









Width of seating area 52 cm | Width of shoulder area (inside) 52 cm

Casada International GmbH
Obermeiers Feld 3
33104 Paderborn | Germany
www.casada.com



# **SPECIALS**



# **BRAINTRONICS®**

Braintronics technology synchronises your brainwaves through audio stimulation. Combined with a specially developed massage programme your stress-ridden body will be able to relax. No matter whether you'd like to relax, dream or learn – anything's possible.




Shiatsu (finger press) is a form of body therapy originating in Japan. The tations. The aim of Shiatsu is to simulate individual parts of the body and to



In this new technique, the massage robot not only moves upwards and downwards, but also forwards and backwards.



Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.




## KNEAD & TAP (DUAL ACTION)


The techniques of kneading and tapping are combined in order to relieve tension and stimulate blood circulation at the same time.



Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles. CARBON HEATING



flow increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain SOLE ROLLER



AIR COMPRESSION Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.



Thanks to the "SL" curve rail technology, which has been approved from human engineering, the range of the massage technology is up to 135cm. This method allows you to enjoy a personal full-body massage from head to thigh and supports the maintenance of the "S"-shaped spine.



For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!

You put yourself in the position of an astronaut – that "weightless feeling".




Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd like.



### BACK AREA An individual applicable to any size massage area.

ADJUSTABLE FOOTREST The footrest can be extended, this way it is adaptable to different body sizes.



### RECLINING POSITION The adjustability of the backrest and footrest contributes to personal well-be-

ADJUST SHOULDER The massage heads can be moved upward or downward to a suitable posi-



# AIR INTENSITY

Adjustable airbag massage intensity.



Whole body airbag massage. Different areas are individually selectable.



The massage heads can massage a particular spot.



Massage speed is adjustable.



The distance between the both massage heads can be set (narrow, medium



## Further relaxation and fun with the music function.



# Wireless data transmission.



### HUMAN HANDS FEELING Specially designed massage heads imitate the hand motion of a professional

masseur. The maximum stretch space is 6,5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebra thoracic vertebra and lumbar vertebra.



While you relax, the integrated USB port can supply, for example, your phone with power.



QUICK START

USB-CHARGER

With one touch the quick start buttons in your armrest allow you to control the most important functions



Arthrosis Traction consists of gentle stretching pressure. It has a stimulating effect on the cartilages, which enhances the transport of bodily fluids, nutrients absorbation and can alleviate symptoms of arthrosis.



Airbags in the lower area of the back ensure gentle stretching of the pelvis. Effectively counteracts tensions.



You can expect a profound full body massage from the neck down to the buttocks.



Using the intelligent memory function the massage chair remembers your favourite massage exactly and begins immediately without a scan.



### ADJUSTABLE FOOTREST The footrest can be extended, this way it is adaptable to different body sizes.

If necessary, the double cushion can be placed to reduce the neck and shoulder massage intensity.



Targeted stretch massages activate the body and have an invigorating effect on the whole body.



### Targeted stretch massages activate the body and have an invigorating effect on the whole of the body



### **ACUPRESSURE POINTS** Key acupressure points are massaged.

all organs inside your body.

Reflexology regenerates foot muscles and has a harmonizing effect on almost



Polluted air around you will be neutralized and you will breathe purified oxygen during your massage.



brainwave stimulation system

The human brain is a complex organ consisting of millions of little nerve cells, called neurons. Each cell constantly produces tiny electrical impulses. Every thought, every sound and every impression, everything is transmitted and communicated by those small electrical impulses. Electricity is, so to say, the language of the brain. The sum of those impulses can be detected and diagrammed, e.g. by using an EEG. Charted on a graph this electrical activity forms a wave-like pattern, know as brainwaves.

With braintronics®, combining audio stimulation and guided meditations with a soothing massage, these brainwaves can be stimulated in order to reach different states of relaxation in the brain. Whether you want to sustainably enhance your concentration, or just relax and pamper yourself - everything is possible with a massage chair from the Casada family.



CASAUL



