

Complete RevitaLens Multi-Purpose Disinfecting Solution

INGREDIENTS

Alexidine 0.00016%; Polyquaternium-1 0.0003%; Boric Acid; Sodium Borate Decahydrate; Tetronic 904; Edetate Disodium (EDTA); Trisodium Citrate Dehydrate; Sodium Chloride; Purified Water.

ALLERGY INFORMATION

N/A

DIETARY

N/A

STORAGE

Store in a cool, dry place. Keep out of sunlight.

MANUFACTURER & ADDRESS

Block B, Liffey Valley Office Campus, Quarryvale, Co Dublin, Ireland.

COUNTRY OF ORIGIN

UK

DISCLAIMER

Please note that whilst we take every care to ensure the product information displayed on our website is correct, product recipes are regularly changed and this may affect nutrition and allergen information therefore you should always check product labels and not rely solely on the information presented here.